

APESIN spray F



Alcohol-based rapid disinfection, liquid - perfume-free

- Fast-acting
- Free of dyes and perfume

Product Profile

- APESIN spray F is a fast-acting disinfectant against bacteria, yeasts, and viruses (limited spectrum).
- This product fulfils the latest hygiene requirements and complies with EU Biocidal Products Regulation (BPR).
- It is suitable for all alcohol resistant surfaces and furniture in healthcare facilities, in industrial and institutional areas as well as in food areas.
- APESIN spray F saves one working step because food contact surfaces do not need to be rinsed with water after the use of the product.
- APESIN spray F is free of chlorine, aldehyde, and phenol.

Area of application

- The product can be used on mattresses, beds, transport trolleys, instrument deposits, carpets and bathroom tiles, dustbins, toilet seats, frameworks, door handles, furnitures and work surfaces, handrails, shoes, as well as in places that are difficult to access.
- Also suitable for the cleaning of cutter boxes and containers. It is best appropriated for knives, boots and cleaning machines.
- APESIN spray F can be applied on aluminium surfaces.
- Do not apply on alcohol sensitive materials (e.g.: PLEXIGLAS®).
- For painted surfaces, test material resistance on an inconspicuous area first.

Application note

- Apply undiluted. Spray directly onto objects from a distance of 20-30 cm, so that surfaces are completely wetted, or apply on a clean and absorbent cloth and wipe the surface completely wet.
- Max. 50 ml/m². Do not spray larger surfaces (max. 2m²).
- Allow the solution to work, for detailed information and exposure time, see next page: Standards & spectrum of activity.

Ingredients:

Ingredients per 100g: 40g ISOPROPYL ALCOHOL, 25g PROPYL ALCOHOL

Sustainable development and production:



Product safety, storage and environmental instructions

Safety: For professional use only. For further information, see Safety Data Sheet. Use biocides safely. Always read the label and product information before use.

Storage: Store at room temperature in the original container.

Environmental: Only give completely emptied container to special waste disposer.

Application & Dosing



Follow below instructions.



Conc.

For disinfection:
Suitable for areas with specific hygiene requirements.



Conc.

Spray application:
Spray onto a clean cloth and wipe the surface wet or spray only.



Surface cleaning:
Wipe the surface wet or damp with a clean cloth.



Conc.

Toilet cleaning:
Clean toilet seats.



Conc.

Ideal for application on surfaces in the kitchen area.




Sales units:

Order No. 714867 10 x 750 ml
Order No. 716058 5L

Your trusted partner on site:



APESIN spray F



Alcohol-based rapid disinfection, liquid - perfume-free

Surface Disinfection: Food area / industrial & institutional area	contact time			
	30 sec.	1 min.	2 min.	5 min.

bactericidal activity				
EN 1276	high organic load		X	
EN 13697	high organic load		X	
yeasticidal activity				
EN 1650	high organic load		X	
EN 13697	high organic load		X	

Surface Disinfection: Medical area	contact time			
	30 sec.	1 min.	2 min.	5 min.

bactericidal activity				
EN 13727	high organic load	X		
EN 13697	high organic load		X	
EN 16615	high organic load		X	
tuberculocidal activity				
EN 14348	high organic load	X		
yeasticidal activity				
EN 13624	high organic load	X		
EN 13697	high organic load		X	
EN 16615	high organic load		X	


virucidal activity				
EN 14476 Virucidal activity against enveloped viruses	high organic load	X		
EN 14476 Limited spectrum virucidal activity	high organic load			X
EN 14476 individual test organisms				
Vaccinia Virus	high organic load	X		
Murine Norovirus (MNV)	high organic load		X	
Adenovirus	high organic load			X

VAH-Certificate	contact time			
	30 sec.	1 min.	2 min.	5 min.

Surface disinfection in medical areas, bacterial and yeasticidal activity	with mechanical action dirty conditions			X
---	---	--	--	---



APESIN spray F



Alcohol-based rapid disinfection, liquid - perfume-free

Surface Disinfection:		contact time			
Food area / industrial & institutional area		30 sec.	1 min.	2 min.	5 min.
bactericidal + yeasticidal activity					
EN 1276 + EN 1650 + EN 13697	high organic load		X		
Surface Disinfection:		contact time			
Medical area		30 sec.	1 min.	2 min.	5 min.
bactericidal + yeasticidal activity (without mechanical action)					
EN 13727 + EN 13624+ EN 13697	high organic load		X		
bactericidal + yeasticidal activity (without mechanical action)					
EN 13727 + EN 13624 + EN 16615	high organic load		X		
tuberculocidal activity					
EN 14348	high organic load	X			
virucidal activity					
EN 14476 Virucidal activity against enveloped viruses	high organic load	X			
EN 14476 Limited spectrum virucidal activity	high organic load				X
EN 14476 individual test organisms					
Vaccinia Virus	high organic load	X			
Murine Norovirus (MNV)	high organic load			X	
Adenovirus	high organic load				X
VAH-Certificate		contact time			
		30 sec.	1 min.	2 min.	5 min.
Surface disinfection in medical areas, bacterial and yeasticidal activity	with mechanical action dirty conditions				X

