



Lobby dust-pan with lid and wheels - White

Hygienic lobby dust-pan for a safe and efficient waste management. The ergonomic design and the functional accessories with it's equipped with, aluminium handle Ø 22 mm, back wheels, lid and bristles-cleaning comb, make it a professional solution in the food sector.



| Code | Colour | Size (mm) | Material | Net weight (kg) |
|--------------|--------|-----------------|----------|-----------------|
| 1082W | White | 270 x 330 x 970 | PP - AL | 0,8 |

Usage information

| | |
|---------------------------------|--|
| Minimum and maximum temperature | min -20 / max 100 C° |
| Washing (dishwasher) | - |
| Autoclave | 121° C for 20 min |
| User and maintenance manual | <ul style="list-style-type: none"> - Remove packaging and any label from tools and disinfect them. - Both manual and autoclave disinfection is possible. - Sanitize tools using appropriate chemical decontamination products, in the concentrations, times and temperatures provided. - Make sure tools completely dry up before next use. - Hang tools while drying up to avoid any contamination or breakage. - The storing area must be clean and far away from sources of intense light or heat. - Do not use damaged tools. |

Compliance

- REACH
- FDA
- EU 1935/2004 and updates
- EU 2023/2006 and updates
- EU 10/2011 and updates

Logistical information

| | | | |
|-------------------------|-----------------|------------------------|----------------|
| Internal packaging type | - | Product barcode | 8008805108244 |
| External packaging type | Carton | Box barcode | 08008805108244 |
| Quantity per pack (pcs) | 6 | Customs tariff | 96039099 |
| Box gross weight (kg) | 5,6 | Country of origin | Italy |
| Box dimensions (mm) | 820 x 355 x 355 | Palletising (pcs) | 90 |
| Box volume (m³) | 0,103 | Pallet dimensions (cm) | 120 x 82 x 192 |

Available versions

| | | | |
|--------|---|----------|--|
| Size | | Material | |
| Colour | B R Y G O | | |